



SR-100 Instructions for Use

1. Go to the App store. Touch the Search Icon. Type in Accelerera, touch Search. You should see Accelerera-Control your Accelerera Device. Download the app. You may need your Apple ID password to load this app.
2. Open the box and locate the wireless charging unit. Locate the charging cable (under the box insert) and plug it into the charging unit and power supply. Plug the power supply into a standard outlet.
3. Locate the Accelerera controller. Make sure the On/Off switch is in the On position.
4. Place the controller on the wireless charging unit. A blue light will come on indicating the unit is charging. This should take about 1.5-2 hours.
5. To use the Accelerera SR-100, remove the controller from the charger and place it in the pouch on the wrap. Connect the actuator lead to the magnetic connector on the controller and place the On/Off switch in the ON position.
6. Position the wrap on your body as instructed by your clinician.
7. Launch the Accelerera app on your phone.
8. You should see a device listed. You can rename the device as you wish, i.e. Left wrist, right knee, personal name.
9. Touch the device name and the control screen should appear. The Stimulation level is controlled by the slider button. Slide it to the right until you can barely feel the vibration. If there is no vibration, make sure there is a blue light lit on the magnetic connector. If not, reposition the magnet to connect to the controller again and use the slider again.
10. Touch Begin. The timer will start documenting the length of use.
11. If you need to pause for any reason, simply touch Pause and then touch Resume when you are ready to continue.
12. When you are finished with your session, place the power unit on the charger.